

Autoclave, Class B, Benchtop Type, STB-BZ Series



1. Automatic door switch
2. Chamber volume: **18L, 23L, 45L**
3. Design temperature: **139°C**
4. Design pressure: **0.25 Mpa**
5. 3 pulsating vacuum, Multiple test function
6. Built-in steam generator

Description

- Automatic door switch
- Quality stainless steel inner tank, sealed cover
- Microcomputer control
- Safety interlocking device
- With BD test, HELIX test, Vacuum test.
- LED display and inductive operation
- Equipped with standard test interface
- Open-type top water tank
- European CLASS B standard.
- With automatic protection function: over temperature protection, over pressure protection, low water level protection, anti dry burning, safety door lock.
- Three times per-vacuum
- Buzzer reminder after sterilization, automatic stop.
- Automatic discharge of cold *air*, automatic steam exhaust after sterilization.
- Built-in steam generator to provide saturated steam quickly
- Built-in printer

Specifications

| Model | STB-B18Z | STB-B23Z | STB-B45Z |
|---------------------------|---------------------|-----------------|-----------------|
| Volume | 18L | 23L | 45L |
| Power | 2.2kw | 2.2kw | 3.5kw |
| Voltage | 110V/220V,60HZ/50HZ | | |
| Design pressure | 0.25Mpa | | |
| Design temperature | 139°C | | |
| Rated working pressure | 0.22Mpa | | |
| Rated working temperature | 134°C | | |
| Chamber volume | Ø250×370 | Ø250×470 | Ø316×615 |
| Outer dimension | 505×620×450 | 505×710×450 | 657×870×575 |