

Class B Autoclave, STB-BB Series



1. Table top Class B Autoclave
2. Chamber volume: 8L, 12L, 18L, 23L
3. Sterilization temperature: 134°C
4. Design pressure: 0.25 Mpa
5. 3 pulsating vacuum, Multiple test function
6. Built-in steam generator

Description

Class B Autoclave is a high level steam sterilizer that uses a vacuum pump to remove air/steam mixtures from the chamber prior to sterilization cycle beginning. It enables the sterilization of both wrapped and unwrapped instruments, no matter their type and complexity, The Class B tabletop autoclave has been the most advanced “plug-n-play” model on the market today. STB-BB Series Class B Autoclave is able to efficiently and effectively sterilizer any load type used in the dental, medical, podiatry, chiropodist, tattoo and veterinary markets.

- Quality stainless steel inner tank, sealed cover
- Microcomputer control
- Safety interlocking device
- With BD test, HELIX test, VACUUM test.



- LCD display and touch buttons.
- With standard test interface
- Open-type top water tank
- European CLASS B standard.
- With automatic protection function: over temperature protection, over pressure protection, low water
- Eevel protection, anti dry burning, and safety door lock.
- Three times pre-vacuum
- Buzzer reminder after sterilization, automatic stop.
- Automatic discharge of cold air, automatic steam exhaust after sterilization.
- Built-in steam generator to provide saturated steam quickly

Specifications

STB-BB Series Class B Autoclave

Model	STB-B8B	STB-B12B	STB-B18B	STB-B23B
Volume	8L	12L	18L	23L
Power	1.5kw	1.8kw	2.2kw	2.2kw
Voltage	110V/220V,60HZ/50HZ			
Design pressure	0.25Mpa			
Design temperature	139°C			
Rated working pressure	0.22 Mpa			
Temp accuracy	±1°C			
Rated working temperature	134°C			
Chamber volume	Ø170×320	Ø200×355	Ø250×370	Ø250×452
Outer dimension	445×580×400	445×580×400	480×640×450	480×690×450